


CPM PACKING INSTRUCTIONS

- REMOVE ALL PADS BEFORE BOXING CPM -

FLIP THE FOOTPLATE ON THE OPTIMA



CORD PLACEMENT FOR OPTIFLEX



WRAP UP THE CORDS



WRAP UP THE CORDS



SURGI-CARE INC. 888-291-8901

CPM PACKING INSTRUCTIONS

- REMOVE ALL PADS BEFORE BOXING CPM -

PLACE THE CPM INSIDE THE BAG PROVIDED, AND SEAL WITH THE RED TAG AND ZIP TIE



PLACE THE CPM IN THE BOX, AND USE THE PACKING MATERIAL PROVIDED TO SECURE CPM



USE THE TAPE ROLL PROVIDED TO CLOSE THE BOX, AND ATTACH RETURN LABEL TO BOX.

**CALL SURGI-CARE 888-291-8901 TO
ARRANGE FOR A PICKUP.**